



[All Collections](#) > [DoseMeRx's Key Features Explained](#) > [Continuous Infusion and DoseMeRx](#)

Continuous Infusion and DoseMeRx

Simulate continuous infusion in DoseMeRx



Written by Sharmeen Roy, PharmD, BCPS

Updated over a week ago

Simulating the outcome of continuous infusion for your patient is quick and easy in DoseMeRx.

In the customize section, simply change the infusion length and the dosing interval to be the same value. For example, an infusion length of 24 hours and a dosing interval of 24 hours would produce a recommendation suitable for continuous infusion.

✓ Customize

Target

Dose

View Outcome

DoseMeRx can calculate a dose based upon previously-processed pathology results. Please review against the guideline dose and consider the clinical implications carefully.

Target

Trough Only

Target Trough

15

mcg/mL

Infusion Length

24

hours

Dosing Interval

24

hours

Number of Doses

6

doses

Next Dose At

2018-09-10 18:00

1

Note: The next dose can only be created from Sep 9 2018, 20:00 to Sep 12 2018, 23:59.

Calculate

Custom Individualized Target Dose Profile [Target is Trough: 15 mcg/mL]

The plot above is generated from the literature model. The recorded pathology results have been processed to calculate a dose-response relationship, which is shown above.

Did this answer your question?

doseme-rx.com

<https://help.doseme-rx.com/en/articles/2316801-continuous-infusion-and-dosemerx>

2/2