



All Collections > Get started with DoseMeRx > Customize a dose or target in DoseMeRx

Customize a dose or target in DoseMeRx

Change your desired target or dose using DoseMeRx's customize feature



Written by Sharmeen Roy, PharmD, BCPS

Updated over a week ago

One of the many popular features of DoseMeRx is to customize your dose or target to determine the best dosing regimen for your patient. You can easily customize the dose, or the target in seconds after [calculating the next dose](#).

Simply navigate to the panel below the dose recommendations to view the customize feature. Enter your custom target or dose parameters and simulate the outcome.

Customize



Target

Dose

DoseMeRx can calculate a dose based upon previously-processed pathology results. Please review against the guideline dose and consider the clinical implications carefully.

Target

Trough Only ▾

Target Trough

15

mcg/mL

Infusion Length

1.5

hours

Dosing Interval

12

hours

Number of Doses

6

doses

Next Dose At

2018-09-10 18:00



Note: The next dose can only be created from Sep 9 2018, 20:00 to Sep 12 2018, 23:59.

Calculate

Tip: If you're logged into DoseMeRx select Help from the top menu, where you can watch our onboarding videos which quickly show you how to customise a dose or target.

Did this answer your question?



doseme-rx.com

