



[All Collections](#) > [Get started with DoseMeRx](#) > [Adding and viewing patient courses](#)

Adding and viewing patient courses

How do I add a course and view my patient's courses?



Written by Sharmeen Roy, PharmD, BCPS

Updated over a week ago

After entering your patient information you have the ability to add a course to your patient from the dropdown menu.



← Add a patient

Patient ID

1234567890

Last Name

Washington

First Name(s)

George

Date of Birth

09 / 25 / 1967

Sex

M

F

Weight

80

kg

Height (OPTIONAL)

176

cm

Unit

Not Recorded

Add a course

Drug Model (OPTIONAL)

Select Drug Model

[Save Patient](#)

If you have selected to add a course with the patient, saving the patient will direct you to the historical plot tab for your patient. If you have not selected to add a course, saving the patient will direct you to the patient course page.

To view your patient courses, simply select a patient to open the patient preview then select 'View all courses'.



PATIENTS

DRUG INFORMATION

ADMIN

HELP

Patients

Patients at Training ▾

+ Add patient

Active ○

Archived □

Search patients



Add filter +

Last 7 days ✕

ID	Name	DOB	Last Modified
<input type="checkbox"/> 1234567890	George Washington	09/25/1967	Aug 19 2020 10:44
<input type="checkbox"/> VANCOMY...	(Pediatric) Two-compartment	01/14/2016	Aug 19 2020 10:02
<input type="checkbox"/> VANCOMY...	(Adult) Two-compartment	12/16/1980	Aug 19 2020 10:02
<input type="checkbox"/> Demo5	Adult Warfarin	06/21/1970	Aug 19 2020 10:02
<input type="checkbox"/> Demo1	Adult Vancomycin	06/17/1957	Aug 19 2020 10:02
<input type="checkbox"/> VANCOMY...	(Neonate) Term	08/02/2020	Aug 18 2020 08:30
<input type="checkbox"/> VANCOMY...	(Neonate) Preterm	08/14/2020	Aug 18 2020 08:30
<input type="checkbox"/> VANCOMY...	(Adult) Hemodialysis	10/04/1954	Aug 18 2020 08:30
<input type="checkbox"/> VANCOMY...	(Adult) Obese	03/16/1984	Aug 18 2020 08:30
<input type="checkbox"/> VANCOMY...	(Adult) One-compartment	12/16/1980	Aug 18 2020 08:30

Patient Preview

George Washington ●

Edit

Archive

ID: 1234567890

DOB: 09/25/1967

Sex: M

Height: 176cm

Weight: 80kg

Unit: Not Recorded

Most Recent Courses: + Add course

Vancomycin IV - Adult (1-comp.)

Add Result

Add Doses

Tobramycin IV - Adult & Pediatric

Add Result

Add Doses

Amikacin IV - Adult

Add Result

Add Doses

[View all courses \(3\) →](#)

Clinicians: Mr. Michael Ferrelli

Archive Multiple Patients

< 1 2 3 >

The courses page will display your patient information at the top of the page. To add a new course simply select, 'Add course' above the course list.



PATIENTS

DRUG INFORMATION

ADMIN

HELP

Patients > George Washington (1234567890)

ACTIVE ●

George Washington ID: 1234567890

Edit

DOB: 09/25/1967 (52 years)

Sex: M Height: 176cm Weight: 80kg

Unit: Not Recorded

Clinicians: Mr. Michael Ferrelli

Created: Aug 19 2020 10:44

Last Updated: Aug 19 2020 10:44

Hide

Courses



Add course

Drug Model	First Dose	Last Dose
<input type="checkbox"/> Vancomycin IV Adult (1-compartment)	No doses recorded.	
<input type="checkbox"/> Tobramycin IV Adult & Pediatric	No doses recorded.	
<input type="checkbox"/> Amikacin IV Adult	No doses recorded.	

Delete Multiple Courses

< 1 >

Course Preview

Vancomycin IV

Adult (1-compartment)

View Course

Delete

Recent Doses: + Add Doses

Dose Time	Interval	Amount
No doses recorded		

[View all doses](#)

Recent Lab Results: + Add Lab Result

Measurement Time	Type	Level
No lab results recorded		

[View all lab results](#)

Did this answer your question?





doseme-rx.com

