

Search for articles...

All Collections > DoseMeRx's Key Features Explained > Reviewing alternative doses

Reviewing alternative doses

Access the feature in DoseMeRx with alternative lower and higher doses to the Individualized or Population recommendations.




Written by Sharmeen Roy, PharmD, BCPS Updated over a week ago


DoseMeRx has a helpful feature that will show you a dose rounding size (e.g. 250mg) higher and lower than the 'Default Individualized Target' recommendation along with the predicted outcomes.

You can use this feature for any of our vancomycin models, once you have <u>added your patient</u>, doses and lab results into DoseMeRx and are ready to <u>calculate the next dose</u>,

To see this feature on the dose recommendation page, simply click on the preview icon.



Here you will see the alternate dose amounts (one size higher and one size lower) over the same dosing interval and number of days.



Simulating one of these alternate doses is easy, scroll down to the '<u>Customize' section</u>, input the new dose amount and simply click 'Calculate.'



Did this answer your question?









doseme-rx.com





